3s PowerCure

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





Quelle: U. Blunck, Marginal adaptation study, *Study Report*, Berlin, 2018.

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
Quelle: C. Arrais, Analysis of pulp temperature and inflammatory response to radiant exposure from an experimental Polywave® LED light curing unit, *Research Report*, State University of Ponta Grossa, 2019.

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Dr. Nathaniel Lawson
Assistant Professor and Division Director of Biomaterials
UAB School of Dentistry, Birmingham, USA

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Quellen:

- [1] N. Ilie, Characteristics of composite and curing unit, *Study Report*, Munich, 2019.
- [2] N. Ilie, Characteristics of composite and curing unit, *Presentation*, Munich, 2018.
- [3] W. Palin, Polymerization characteristics of Tetric EvoCeram Bulk Fill and F-Composite 2, *Study Report*, Birmingham (UK), 2015.
- [4] W. Palin, M. Hadis, High irradiance polymerization of "flash-cured" resin composites, *Study Report*, Birmingham (UK), 2018.
- [5] S. Lenz, Research Report of Tetric PowerFill, *Test Report*, Ivoclar Vivadent, 2019.
- [6] B. Gebhardt, Verification Report Tetric PowerFill, *Test Report*, Ivoclar Vivadent, 2017.
- [7] S. Heintze, Design validation report Tetric PowerFill, *Test Report*, Ivoclar Vivadent, 2017.
- [8] R. Hirata, Effect of sonic resin composite delivery on void formation assessed by micro-computed tomography, *Operative Dentistry* 2018, 43-2, 144-150.
- [9] M. Cowen, J. M. Powers, Bond strength laboratory evaluation of a new curing light and restoratives, *Study Report*, Dental Advisor, 2018.
- [10] T. Bock, Research report Adhese Universal, Shear Bond Strength (SBS), *Research Report*, Ivoclar Vivadent, 2018.
- [11] U. Blunck, Marginal adaptation study, *Study Report*, Berlin, 2018.
- [12] A. Lebedenko, Comparative fillings: Conventional layering technique versus 3sCure two-layer technique, *Test Report*, Ivoclar Vivadent, 2018.
- [13] L. Enggist, Comparison of conventional workflow in direct filling therapy with 3sCure System: Time need in three different German dental offices, *Test Report*, Ivoclar Vivadent, 2020.
- [14] N. Ragazzini, Comparison of restorative time of direct fillings class I & II placed with traditional layering technique or bulk layering technique, *Clinical Report*, Bologna, 2020.
- [15] N. Lawson, Clinical evaluation of a bulk fill resin composite, 24-month report, *Study Report*, Birmingham (US), 2020.
- [16] B. Ganster et al, New photocleavable structures, 4. Acylgermane-based photoinitiator for visible lightcuring, *Macromolecular Rapid Commun.* 2008, 29, 57-62.
- [17] B. Ganster et al, New photocleavable structures, Diacylgermane-based photoinitiator for visible lightcuring, *Macromolecules* 2008, 41, 2394-2400.
- [18] S.E. Kopperud et al, Light curing procedures – performance, knowledge level and safety awareness among dentists, *J. Dent.* 2017, 58, 67-73.
- [19] F. C. Calheiros, et al, Degree of conversion and mechanical properties of a BisGMA:TEGDMA composite as a function of the applied radiant exposure. *J Biomed Mater Res B Appl Biomater* 2008, 84, 503-509.
- [20] C. Arrais, Analysis of pulp temperature and inflammatory response to radiant exposure from an experimental Polywave[®] LED light curing unit, *Research Report*, State University of Ponta Grossa, 2019.
- [21] Bolis C et al., *Oral Health Prev. Dent.* 2015, 13, p. 545-556.

Hersteller und Vertrieb
Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclar.com

Vertrieb Deutschland
Ivoclar Vivadent GmbH
Dr. Adolf-Schneider-Str. 2
D-73479 Ellwangen, Jagst
Tel. +49 7961 8890
Fax +49 7961 6326
www.ivoclar.com

704397/DE/2022-07-04

ivoclar
vivadent®