Recommended grinding instruments

	Firing		Polishing		
	IPS e.max® CAD	IPS e.max® ZirCAD	IPS Empress® CAD	Tetric® CAD	Telio® CAD
	Lithium disilicate glass-ceramics (LS ₂)	Zirconium oxide ceramics (ZrO ₂)	Leucite glass-ceramics	Composite	PMMA
Smoothing out the attachment point	 Fine diamond burs or silicon carbide burs ^{[1], [5]}	 Fine diamond burs or fine tungsten carbide burs ^{[1], [5]}	 Fine diamond burs or silicon carbide burs ^{[1], [5]}	 Fine diamond burs ^{[1], [5]}	 Cross-cut tungsten carbide burs ^[1]
Finishing	 Fine diamond burs or silicon carbide burs ^{[1], [5]}	 Fine diamond burs or fine tungsten carbide burs or silicon carbide burs ^{[1], [5]}	 Fine diamond burs or silicon carbide burs ^{[1], [5]}	 Fine diamond burs ^{[1], [5]}	 Fine diamond burs or cross-cut tungsten carbide burs ^[1]
Polishing	 2-step polishing using OptraGloss® ^{[1], [3]}	 Optional: 1-step polishing using OptraGloss® ^{[1], [3]}	 2-step polishing using OptraGloss® ^{[1], [3]}	 1-step polishing using OptraGloss® ^{[1], [3]}	 1-step polishing using OptraGloss® ^{[1], [3]}
Cristallization or sintering	 Programat® CS6		General notes <ul style="list-style-type: none"> - The respective indications of the manufacturer of the grinding tools on the proper use, e.g. speed, have to be observed. - Please use light pressure during processing irrespective of the material to avoid overheating. 		
Corrections	 Extra-fine diamond burs ^{[1], [5]}	 Fine diamond burs or silicon carbide burs ^{[1], [5]}	^[1] Dry processing ^[2] Wet processing ^[3] Rotary speed max. 10,000 rpm ^[4] Rotary speed max. 15,000 rpm ^[5] Rotary speed max. 20,000 rpm		
Polishing	 2-step polishing using OptraGloss® ^{[1], [3]}	 2-step polishing using OptraGloss® ^{[1], [3]}			

Recommended grinding instruments

	Firing			Polishing	
	IPS e.max® CAD	IPS e.max® ZirCAD	IPS Empress® CAD	Tetric® CAD	Telio® CAD
	Lithium disilicate glass-ceramics (LS ₂)	Zirconium oxide ceramics (ZrO ₂)	Leucite glass-ceramics	Composite	PMMA
Extensive corrections	Fine diamond burs [2], [4]	Fine diamond burs [2], [4]	Fine diamond burs [2], [4]	Fine diamond burs [2], [4]	Fine diamond burs [2], [4]
Minor corrections	Extra-fine diamond burs [2], [4]	Fine diamond burs [2], [4]	Extra-fine diamond burs [2], [4]	Extra-fine diamond burs [2], [4]	Extra-fine diamond burs [2], [4]
Polishing	2-step polishing using OptraGloss® [2], [3]	2-step polishing using OptraGloss® [2], [3]	2-step polishing using OptraGloss® [2], [3]	1-step polishing using OptraGloss® [2], [3]	1-step polishing using OptraGloss® [2], [3]
Endo access	Medium diamond burs	Medium diamond burs	Medium diamond burs	Medium diamond burs	-

The necessary processing steps depend on the chosen processing method and are to be observed. For further information please see the Instructions for Use of the individual processing methods.

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Ivoclar Vivadent AG