

Versatile applications

SR Ivocron®

High-quality
PMMA veneering material



One material for three techniques

SR Ivocron is based on many years of clinical experience. The high-quality PMMA veneering material is suitable for diverse areas of application. Thanks to the individual choice of processing technique, it can be used to create esthetic composite veneers on crowns and bridges, for temporary applications and for the fixation of denture teeth in prosthetics.

Pouring technique

Cold technique

- ✓ Short- and medium-term temporaries
- ✓ Filling and fixation of ground denture teeth on model cast frameworks
- ✓ Repairs

Layering technique

Hot technique

- ✓ Crown and bridge veneering technique
- ✓ Long-term temporaries

Flasking technique

Press technique

- ✓ Crown and bridge veneering technique
- ✓ Long-term temporaries



- ✓ Esthetic restorations in the 20 shades of the Chromascop shade guide and in Gingiva pink
- ✓ Compatible with PMMA denture teeth
- ✓ Optimized bond to metal frameworks in conjunction with the SR Link bonding agent

Ivoclar Vivadent AG | Bedererstrasse 2 | 9494 Schaan | Liechtenstein
Tel. +423 235 35 35 | Fax +423 235 33 60

[ivoclar.com](https://www.ivoclar.com)